



JANAK[®]
Because we care.



Bx2000 (Manual Bed)

User's Manual

Index



1.	Safety
2.	Definitions
3.	Introduction
4.	Intended Use
5.	Installation
6.	Operation
7.	Optional Accessories
8.	Maintenance & Cleaning
9.	Warranty & Service

SAFETY

In this instruction manual the following safety symbols are used



Danger: This symbol will appear wherever safety instruction is designed to protect people from physical harm. The symbol stands from imminent danger of death or serious injury.



Caution: This symbol will appear where situations are described, which might be dangerous, and which might inflict slight injuries.



Advice: This symbol will appear in front of additional helpful pieces of advice.

- A dot in front of the text means, this is what you have got to do.

Safety Instructions

- Always use accessories that have been designed and approved for use with the bed.
- Ensure that the bed is assembled and installed in accordance with the instructions given in this user manual.
- Do not use this bed for transporting patients.
- Restrict child visitors without proper attendants. Supervise the bed to avoid any mishap due to unintended movement of the bed.
- Do not allow any visitor to sit on the bed, as it may cause overloading. (Refer section - Operation).
- Don't move the vacant bed in a vertical condition or push on the floor, as this will damage the leveling & alignment of the hospital bed. To shift the bed, lift the bed from the main frame on both the sides and then move.
- Accessories or any devices used along with the bed should be as per its intended use specified and approved by its manufacturer.
- Keep the user manual handy with the operator/user for ready reference.
- Always use a mattress of the correct size and type. Before operating the bed, make sure that the patient is safely positioned to avoid entrapment or imbalance.



Use of water mattress is prohibited.

DEFINITIONS

Operator

Operator is the person appointed by the owner (i.e. clinic, hospital management, etc.) who will be responsible for the safe operation of this hospital bed.

User

User is a person who is capable of assessing clinical condition of patient e.g. specialist medical personnel, doctors, nursing staff, attendants etc. On the basis of their training experience or through the instructions, they are entitled to operate the bed on their own responsibility or to carry out work on it and are able to recognize and control possible risks to the patient as well as the bed.

Patient

The patient is a person who is ill, infirm, disabled or in need of care & occupying this hospital bed. Each time the bed is allotted to a new patient, it is recommended that the patient is instructed by the operator or the user about the functions important for his/her comfort and safety.

INTRODUCTION

About this user manual



This user manual contains instructions for the proper and safe use of the bed and for the patient safety and comfort, which has to be strictly followed.

This user manual has to be read and observed by every person using or handling the bed. In case of any questions, kindly contact **Midmark**

INTENDED USE

- These manual beds are non-invasive in nature and shall be used for resting of the patients by medical/paramedical/hospital authorities to help and support them in diagnosis, prevention, monitoring, treatment or alleviation of the disease.
- These manual beds are patient support and patient rest facility, which help medical experts in carrying out their medical procedures comfortably.
- These are manually operated without any electrical or electro mechanically operated functions
- The manually operated beds are manufactured and supplied in non-sterile condition for use in normal hospital environment
- Use bed accessories in accordance with the manufacturer's instructions for patient safety and comfort. Use recommended spares and accessories like safety side rails, saline rod and mattress of suitable design manufactured by "Midmark" to eliminate hazards.
- The hospital bed is designed to be dismantled into sections for ease of transportation, but can be easily assembled with a minimum of tools.

Bed Features



Please refer the product catalogue for this type of product. (Bx2000)

Technical Data



Please refer the product catalogue for this type of product. (Bx2000)

Optional Accessories



Please refer the product catalogue for various accessories available and compatible with this product. (Bx2000)

INSTALLATION

1. Unpack all the components carefully.
2. If the packing is drenched by rain, water or other liquids, inform this matter to the operator.
3. Each bed is supplied with the following components as per packing list -

• Main frame with backrest screw mechanism	01 No
• Legs (As per customer order specification)	02 Nos
• Head & Foot Boards (As per customer order specification)	02 Nos
• Polymer Moulded Handle	01 No
• Safety Side Rails (As per customer order)	01 Pair
• Mattress (As per customer order)	01 No
4. Fix the legs on the main frame.
5. Place Head & Foot Boards on the main frame.
6. Fix the Safety Side Rails in position.
7. Fix the optional accessories on the bed.
8. Demonstrate the functions of the bed to the user/operator.



Installation of the hospital bed should be done only by Midmark trained technicians.



Packing materials has to be disposed off according to the local regulations for waste disposal.



With due consideration to environmental limitations, materials used to make the device are not affected by contact with any bio chemicals or bio waste. When used as per its intended purpose, all materials used in the product and accessories do not facilitate breeding or growth of bio contaminants on their own.

OPERATION

- Before putting the hospital bed into use for the first time, ensure all the transport securing devices and packing film are removed.
- Clean and disinfect the hospital bed prior to the first time use.
- The hospital bed should be used on a flat floor.
- Check, if all the adjustments are functioning properly.
- This hospital bed may only be operated by person who has the knowledge of its operation.



Do not overload the bed above 135 kgs, which is the safe working load of this bed (Patient + Mattress).

$$\frac{\text{Diagram of bed with mattress and handles}}{\text{Warning icon}} = 100 \text{ kg} \quad \frac{\text{Diagram of bed with patient and mattress}}{\text{Warning icon}} = 135 \text{ kg}$$

Backrest



For adjusting the backrest section, first insert the handle (1) on the screw head (2) below the frame on the leg end of the bed, then rotate the handle anticlockwise to raise or clockwise to lower the backrest section (see Figure 1).



Do not turn the handle beyond its limits in any direction, as it may cause damage to the mechanism. Always remove the handle, when not in use or while transporting the bed. It may hurt people passing by or the handle can also be damaged.

Head & Foot Boards (Removable Type):



(0181-Conventional Bow)



(0182-S.S. Laminated Bow)



(0183-Polymer Moulded Bow)

Head Boards:

- The head board is located at the head end of the bed frame.
- The head board can be removed or attached in a single step by extending equal force with both the hands in the bracket by ensuring the proper matching of the brackets with board holder fitted on the main frame.
- The head board can be removed for the ease of access to the patient from head side.

Foot Boards:


- The foot board is located at the foot end of the bed frame.
- The foot board can be removed or attached in a single step by extending the equal force with both the hands in the bracket by ensuring the proper matching of the brackets with board holder fitted on the main frame.




Do not hang any heavy items on the Head & Foot Board.

Side Rail Options


Safety Side Rails (0354, 0356A, 0358)



(0354-Swing Down Type)



(0356A-Collapsible Type)



(0358-Collapsible Type with S.S. Tubes)

0354 - Swing Down Type Safety Side Rails

1. To lock the railing, pull the railing upwards by extending equal force with both the hands. On pulling the railing completely upwards the Lever (1) will lock the railing.
2. To unlock the railing, pull the Lever (1) with one hand & hold the railing with the other hand. On releasing the Lever, place the railing downward.

0356A - Collapsible Type Safety Side Rails (M.S. Powder Coated)

1. To lock the railing, pull up the railing with the hand till the plunger Lock (1) snaps into place in the locking slot.
2. To unlock the railing, pull out the plunger Lock (1) with one hand & hold the railing with the other hand. Place the railing downward completely till it rests on the bed frame.

0358 - Collapsible Type Safety Side Rails (with S.S. Tubes)

1. To lock the railing, pull up the railing with the hand till the plunger Lock (1) snaps into place in the locking slot.
2. To unlock the railing, pull out the plunger Lock (1) with one hand & hold the railing with the other hand. Place the railing downward completely till it rests on the bed frame.

- Safety Side Rails should be properly fitted and secured by a trained technician.
- Always use the Safety Side Rails with the mattress of the correct size, which is approved for use with this bed for the safety of the patient.
- Do not use the Safety Side Rails to move or lift the bed.



Although the Safety Side Railings have been designed to reduce the risk of patient injury, the potential exists for patient entrapment, particularly in agitated or disoriented patients, patients who have cognitive disability or have uncontrolled body movement as well as patients who lack the physical strength to extract themselves, if they get entrapped.

Users should carefully evaluate the need for safety side rails and periodically check patients in accordance with facility protocols for safe positioning.



If these safety measures are not observed, injuries can be caused to hands, knees, fingers, feet, legs or other body parts of patient due to unintended and unintentional entrapment. In children or weak patients who are less than 146 cms tall, non observance of these guidelines may even lead to death.

OPTIONAL ACCESSORIES




IV Pole:

IV Poles are used to hang IV bags for infusion of liquid substances directly into a vein. The IV pole is removable, telescopic pole that can be installed in any of the sockets located at four locations provided on the main frame of the bed. The IV pole has a safe working load of 5 kgs.



Do not exceed the load capacity of the IV pole. If IV pole is overloaded, personal injury or equipment damage may occur.



Mattress



Patient should not be put on the bed without the right mattress.



Midmark recommend the use of Midmark made mattresses that have been designed & tested specifically for the BX 2000 bed. Also ensure that the dimensions & construction of mattress should minimize the gaps where entrapment could occur, provide sufficient height between the mattress & the top of the side rail to prevent accidental roll- over events, provide appropriate firmness at the edges of the mattress to facilitate safe transfers in & out of the bed & do not interfere with the proper operation of the side rails.



Mattress impermeability could be affected by needle sticks. Users should be instructed to avoid punctures caused by improper use of needle sticks. Failure to do so could result in cross-infection & patient injury. The mattress should be regularly inspected for punctures, rips, tears or other such damages & should be replaced with suitable mattress, if required.



Patients should not be allowed to smoke in bed as it may cause fire hazard, and personal injury or equipment damage could occur. Mattress, sheets & pillows generally do not have flame retardant properties.



Urine Bag Holder

Urine bag holder is used to hang the urine bag. The urine bag can be placed on the hooks provided on the urine bag holder.



S.S. Lifting Rod

The lifting rod is used by the patients who have difficulty in getting up from the bed on their own.



Castors with Rubber buffer



The castors are intended for indoor use only. Outdoor use may cause temporary or permanent damage to the castors. The bed is designed to be moved on smooth stepless floor through 1000 mm wide door way.



To prepare the bed for transport (without patient) –

- All fixed equipments and accessories should be secured properly to the bed before transportation.
- All detachable accessories should be safely removed from the bed before transportation.
- The safety side rails should be in lowered position.
- Unlock the brakes on the castors before moving the bed.
- If the bed is not equipped with the castors, it should be lifted from its main frame & transported from one place to another.

MAINTENANCE & CLEANING

Each time before putting the bed into service the user must ensure that -

- All the functions are in satisfactory working condition before allotting the bed to the patient.
- Lubrication of all the mechanical joints/pins/moving parts/tightening of the bolts should be carried out at least once in three months, to avoid unnecessary sound and wear & tear of the bed.
- Visual inspection is carried out at regular intervals and check for any mechanical damage (e.g. loose /missing rivets, screws, etc.)
- If any damage or malfunction is observed, transfer the patient to other bed and report this matter to the operator.
- If the bed is not being used, it should be stored in a dry and flat indoor environment which is rodent free with gradient not more than 5 degrees.
- If the bed is not in use for long time, it is necessary to conduct a test to check for functioning and any mechanical damage that might have been caused.
- Depending on the usage certain moving parts may require replacement due to wear & tear. We recommend timely replacement with genuine 'Midmark' spares for safety and trouble free operation.

Cleaning & Disinfecting -

- The bed should be properly cleaned and disinfected to avoid infection and stains before the bed is allotted to a new patient.
- Immobilize the bed & remove all the detachable accessories and components. Use suitable rubber or plastic sheet below the mattress to prevent dampness from coming in contact with the bed's support surface. Any prolonged contact of urine or acid on mild steel can corrode the surface. Keep the surface dry at all times.
- Wearing suitable protective clothing, clean all surfaces with a disposable cloth soaked in a neutral detergent and hot water. We recommend use of Bacillocid to clean the bed.
- Start by cleaning the upper sections of the bed & work along all horizontal surfaces. Work methodically towards the lower sections of the bed. Take extra care to clean areas that may trap dust or dirt. All detachable accessories and its component can be cleaned separately.



- Do not expose the bed to excessive moisture as corrosion may cause personal injury or equipment damage.
- Do not use any steam cleaning device on the bed. The bed is not suitable for machine wash. Do not pour water on bed or use high-pressure water jets for cleaning purpose.

Disposal



Product and its defective parts especially elements made of plastic, steel & rubber material should be disposed off in accordance with the statutory rules and regulations prevailing for such disposal of waste. The user is advised to contact 'Midmark' for further assistance.

ROUTINE INSPECTION CHECK LIST : INSPECTION BY THE USER

(Recommended Frequency: Every new allotment but not later than one month)

Checks	Observation	OK	Not OK	Description of faults
Visual check of the mechanical components				
Bed frame	Damage/Deformation			
Screws and rivets are tightened	Missing/Loose			
Performance check of all mechanical components				
Detachable head and footboards accessories	Secure retention Fastening, damage, missing			
Sign	Result of Inspection	Date		

MALFUNCTIONS

PROBLEM	POSSIBLE CAUSES	REMEDY
Backrest not operating	Obstruction/Broken screw mechanisms	Replace the screw mechanism

WARRANTY & SERVICE

Midmark Pvt. Ltd. warrants, solely to the end users of this bed that for a period specified in the warranty certificate, such bed shall perform as per its intended use. The foregoing warranties shall not apply if the product has been subjected to: (a) physical abuse, misuse, abnormal use, or use not consistent with Midmark Pvt. Ltd. Ltd.'s published directions, (b) any modification by anyone other than Midmark Pvt. Ltd. or third parties designated by Midmark Pvt. Ltd. ; or (c) fraud, tampering, negligence or accidents, depending upon the nature of the products.

Limitation of Liability: In no event shall Midmark Pvt. Ltd. be liable or obligated to the end user in any manner for any special, non-compensatory, consequential, indirect, incidental, statutory or punitive damages of any kind, including without limitation, for injury or loss of life, lost profits, lost sales, lost revenue, loss of any software or data, or loss of use of any hardware, software or data, regardless of the form of action, whether in contract, tort, negligence, strict product liability, or otherwise, even if Midmark Pvt. Ltd. has been informed of or is aware of the possibility of any such damages in advance. Midmark Pvt. Ltd.'s total aggregate liability under this agreement and in relation to anything which Janak Healthcare Pvt. Ltd. has done or not done in connection with this agreement shall be limited to the price of the products purchased and paid for by the end user during the period in which such liability arises. The limitations set forth above shall be deemed to apply to the maximum extent permitted by applicable law and notwithstanding the failure of the essential purpose of any limited remedies.

In order to maintain functional safety and right to claim under the warranty, only original 'Midmark' replacements parts must be used. For ordering replacement parts, customer service request and further question, please contact our technical support department.

MANUFACTURER'S ADDRESS

Corporate Office:

Midmark India Pvt. Ltd.

Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070

Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100

www.midmark.in | india@midmark.com

Registered Office:

Midmark India Pvt. Ltd.

Janak House, Opp. Indian Oil Corp. Depot, Sheikh Misry Road, Wadala (E), Mumbai - 400 037.

Tel: +91 22 2412 0171/2413

Toll Free No.: 1800 22 8020, Email: technicalsupport@midmark.com



European Representative Address



MIDMARK Europe SAS
22, Rue de, Saint Denis,
53500 ERNEE,
France



Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City,
LBS Road, Kurla (W), Mumbai - 400 070

Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100

www.midmark.in | india@midmark.com

Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships

