



1803C - Health Check Examination Couch

- Overall approximate dimension: 1895 mm x 560 mm x 830 mm, 75" x 22" x 33" (L x W x H). With BP tray: 945 mm W, 37".
- Mattress platform size: 1895 mm x 560 mm x 64 mm, 75" x 22" x 2.5" (L x W x H).
- Upper section of box size: 1216 mm x 460 mm x 630 mm, 48" x 18" x 25" (L x W x H).
- Three cabinet size: 352 mm x 435 mm x 435 mm, 14" x 17" x 17" (W x H x D).
- B. P. tray size: 127 mm x 365 mm x 20 mm, 5" x 14" x 1" (W x L x H).
- Fixed upholstered top in two pieces of 560 mm, 22" (W).
- Headrest adjustable on gas spring.
- Upper section of box with three drawers.
- Lower section comprises of three cabinets with separate doors and lock on two side doors.
- B.P. apparatus tray provided near headrest.
- Inbuilt sliding step stool.
- Pre-treated and powder coated.
- Supplied in Semi Knocked Down condition (SKD).
- Safe working load of 170 kgs and patient load bearing capacity of 135 kgs.

Optional Accessories

Paper roll holder (0175)

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec