



## 1803M - Health Check Single Top Examination Couch with "U" Cut

- Overall approximate dimension: 1860 mm x 560 mm x 830 mm, 73" x 22" x 33" (L x W x H).
- Mattress platform size: 1830 mm x 560 mm x 64 mm, 72" x 22" x 2.5" (L x W x H).
- Upper section of box size: 1216 mm x 460 mm x 630 mm, 48" x 18" x 25" (L x W x H).
- Three cabinet size: 352 mm x 435 mm x 435 mm, 14" x 17" x 17" (W x H x D).
- Single section upholstered top, with "U" cut.
- A pair of synthetic rubber padded S.S. lithotomy rods with rexine ankle straps.
- M.S. powder coated drawers and cabinet doors and sliding inbuilt step stool.
- Upper section of box with three drawers.
- Lower section comprises of three cabinets with separate doors and lock on two side doors.
- B.P. apparatus tray provided near headrest : 127 mm x 365 mm x 20 mm, 5" x 14" x 1" (W x L x H).
- Pre-treated and epoxy powder coated.
- Supplied in SKD condition.
- Safe working load of 170 kgs and patient load bearing capacity of 135 kgs.

### Optional Accessories

Paper roll holder (0175)