

Jx4000



Jx4000 Motorized Gynaecology Chair



Jx4000 Specifications

Overall Approximate Dimensions	
Overall Length - Jx4000	1660 mm 65" (Inclusive of foldable leg section of 360 mm)
Overall Length - Jx4001	1780 mm 70" (Inclusive of removable leg section of 490 mm)
Minimum Height (electrically operated)	555 mm 22"
Maximum Height (electrically operated)	915 mm 36"
Overall Approximate Width	700 mm 28"
Upholstered Platform Width	580 mm 23"
Upholstery Thickness	65 mm 2.56"
Features (Inbuilt)	
Patient Hand Control	User-friendly handset for height adjustment
Backrest	Gas spring assisted
Backrest Angle	0 - 65°
Pelvic Tilt	Manually controlled
Paper Roll Holder	Integrated
Pair of Lithotomy Crutches with Cover	Incorporated
Plastic Pan	Inbuilt
Hand Grips	Provided
Safe Working Load	135 kgs
Patient Load Bearing Capacity	125Kgs
Battery Backup	Available
Electrical Specifications	
Power In (Mains Power Supply)	Nominal 230 V Ac 90 V Ac to 300 V Ac operating range with battery
Protection Voltage	45 Hz to 55 Hz
Electrical Shock Protection	Class I
Degree of Shock Protection	Type B
Liquid Ingress Protection	IP 54
Rechargeable Batteries	2 x 12 Volt Sealed Lead Acid
Duty Cycle	10% (Two minutes for every eighteen minutes)



Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec



We are interested in growing relationships

Jx4000;02:25.10.17

Jx4000



Jx4000



Welcome Position

Flat Position

Examination Position



Hand Grips



Lithotomy Crutches



Plastic Pan



Manual Control of Pelvic Tilt



Foldable Section (for Jx4000)



Removable Section (for Jx4001)

Comfort

- The pelvic tilt ensures a good angle for examining the pelvic examination.
- The low height allows patients to access the table alone, avoiding back pains to doctors and nursing staff.
- Gas spring assisted backrest helps in better sitting comfort and examination.

Safety

- The low height allows easy access, which is especially useful for pregnant ladies.
- Inbuilt battery backup.
- SMPS control box helps in regulating the electrical power efficiently.

Accessibility

- The height adjustment facility allows the doctor to treat the patient more efficiently.
- This adjustment allows the doctor to keep eye to eye contact, ensuring distressful examination.

Hygiene

- Easy to clean.
- Upholstery can be changed when required.
- Integrated paper roll holder helps in maintaining better hygiene.

For upholstery color options, refer color shade card.