

ICU Solutions



0104A Paediatric ICU Bed

Overall Approximate Dimension

Buffer to Buffer	1500 mm x 822 mm, 59" x 32" (L x W)
Platform Dimension	1451 mm x 772 mm, 57" x 30" (L x W)
Mattress Platform	1340 mm x 655 mm, 53" x 26" (L x W)
Buffer to Handle	1680 mm, 66"
Fixed Height	599 mm, 24"

Features (Inbuilt)

Backrest Angle	Atleast 0° - 70°
Trendelenburg Angle	Atleast 0° -10°
Reverse Trendelenburg Angle	Atleast 0° - 10°
Four Corner Rubber Buffer of 100mm dia	Available
S.S. Tubular, Telescopic, Collapsible Side Rails	Available
125 mm High Grade Synthetic Castors two with brake & two without brake	Available
S.S. Laminated Head & Foot Bows	Available
Telescopic I.V. Pole with 2 hooks (0154)	Available
Removable Polymer moulded handle	Available
Powder Coating	Available
Patient Load Bearing Capacity	Safe working load of 110 kgs

Optional Accessories

Two section mattress with 3" thick PU foam of 32 density covered with PVC (0122A)	Available
Urine Bag Holder (0156)	Available
S.S. Telescopic I.V. Pole with four hooks (0154C)	Available
Chart Holder (0158A)	Available