



0104H Paediatric ICU Bed with height

Overall Approximate Dimension

Buffer to Buffer*	1535 mm x 860 mm, 60" x 34" (L x W)
Platform Dimension*	1450 mm x 770 mm, 57" x 30" (L x W)
Mattress Platform*	1330 mm x 680 mm, 52" x 27" (L x W)
Buffer to Handle	1720 mm, 68"
Height Adjustment**	455 mm to 710 mm, 18" to 28"

*(\pm 10 mm Engineering Variation) **(\pm 5 mm Engineering Variation)

Features (Inbuilt)

Backrest Angle	Atleast 0 ^o - 70 ^o
Trendelenburg Angle	Atleast 0 ^o - 12 ^o
Reverse Trendelenburg Angle	Atleast 0 ^o - 08 ^o
Four Corner Rubber Buffer of 125 mm dia	Available
S.S. Tubular, Telescopic, Collapsible Side Rails	Available
125 mm High Grade Synthetic Castors two with brake & two without brake	Available
S.S. Laminated Head & Foot Bows	Available
Telescopic I.V. Pole with 2 hooks (0154)	Available
Removable Polymer moulded handle	Available
Powder Coating	Available
Safe working load	110 kgs
Patient Load Bearing Capacity	90 kgs

Optional Accessories

Two section mattress with 3" thick PU foam of 32 density covered with PVC (0122A)	Available
S.S. Telescopic I.V. Pole with four hooks (0154C)	Available
Oxygen Cylinder Cage (0166)	Available
Chart Holder (0158A)	Available
Urine Bag Holder (0156)	Available