



9003 Crash Cart

- Overall approximate dimension: 955 mm x 493 mm x 1554 mm, 38" x 19" x 61" (L x W x H)
- S.S. Tubular Frame & S.S. Shelves
- Six multicoloured removable bins and two polystyrene lockable storage units with three drawers having containers of different sizes
- Four 125 mm, 5" diameter castor wheels with high grade synthetic body, two with brake and two without brake
- Complete with corner buffers, powder coated oxygen cylinder holder, catheter holder, cardiac massage board, S.S. I.V. rod
- Supplied in KDC

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships

