

Cx4000A Comfort Bed



This bed combines the full features of hospital bed at very economic rate with adjustment of back, upper leg height and trendelenburg / reverse trendelengburg on four separate crank systems. This bed is very useful to clinics, nursing homes and hospitals. With the trademark 'Midmark India' desing and quality, this bed is what we recommended as a minimum hospital bed for the comfort of your patients and nursing staff.

Technical Data

Cx4000A

Overall Approximate Dimension	
Buffer to Buffer	2140 mm x 970 mm, 84" x 38" (L x W)
Platform Dimension	2090 mm x 920 mm, 82" x 36" (L x W)
Mattress Platform	1980 mm x 865 mm, 78" x 34" (L x W)
Buffer to Handle	2325 mm (L)
Minimum Height	490 mm 19"
Maximum Height	725 mm 29"
Features	
Removable Head & Foot Panels	Available
Raised Backrest	Available
Raised Upper Leg	Available
Calf Elevation	Available
Four Section Metallic Perforated Tops	Available
125 mm High Grade Synthetic Castors, Two with Brake and Two without Brake	Available
Backrest Angle	0 - 65°
Knee Rest Angle	0 - 50°
Trendelenburg Angle	0 - 10°
Reverse Trendelenburg Angle	0 - 7°
Buffers at 4 Corners	Available
I.V. Pole Provision	Available
Four separate screw mechanism for Backrest, Upper Leg, Trendelenburg & Reverse Trendelenburg Tilt and Height Adjustments	Available
Polymer Moulded Handle	Available
All M.S. Parts are Pre-treated and Powder Coated	YES
Patient Load Bearing Capacity	100 kgs
Safe Working Load	135 kgs
Option of Head and Footboards	
Polymer Moulded Head and Footboards (0183) (1 Pair)	Available
S.S. Laminated Head and Footboards (182G) (1 Pair)	Available
Optional Accessories	
Four Section Mattress with 4" Thick PU Foam of 40 Density Covered with PVC (0321) - 1980 mm x 865 mm, 78" x 34" (L x W)	YES
Oxygen Cylinder Cage (0166E)	YES
S.S. Heavy Duty I.V. Pole with Two Hooks (0164)	YES
S.S. Heavy Duty I.V. Pole with Four Hooks (0164C)	YES
Telescopic I.V. Pole with Two Hooks (0154)	YES
Telescopic I.V. Pole with Four Hooks (0154C)	YES
Urine Bag Holder (0156)	YES
Moulded Chart Holder (0158A)	YES
S.S. Lifting Pole (0553)	YES
Orthopaedic Traction Pulley (0187A)	YES
Full Length S.S. Collapsible Safety Side Railings (0352)	YES
Full Length Swing Down Type Side Rails EPC (Pair) (0354)	YES
M.S. Collapsible Type Safety Side Rails Powder Coated (Pair) (0356A)	YES
Collapsible Railings with Stainless Steel Tubes (0358) (Pair)	YES