

# Cx4000 PLUS Comfort Bed



Presenting the all new Cx4000 PLUS, a thoughtfully crafted variant designed for the challenging healthcare environments, keeping caregiver's needs in mind!

It is packed with the full set of clinical positions like Height, Backrest, Trendelenburg, Reverse Trendelenburg, Knee Rest on individual smooth crank systems and a multi position DVT end platform, making it a perfect choice for critical care and high dependency areas. Armed with an elegantly designed and smoothly operated split type polymer molded safety side rails, it not only offers safety and hygiene but also brings treatment areas to life. The bed is engineered for trouble free service life with a tough frame design, strong support linkages and high grade synthetic castors.

# Cx4000 PLUS Comfort Bed



## Technical Data

Cx4000 PLUS

Overall Approximate Dimension	
Buffer to Buffer (L x W)	2140 mm x 1010 mm, 84" x 40", ± 10 mm
Platform Dimension (L x W)	2080 mm x 950 mm, 82" x 37", ± 10 mm
Mattress Platform (L x W)	1905 mm x 865 mm, 78" x 34", ± 10 mm
Buffer to Handle (L)	2340 mm, ± 10 mm
Minimum Height	520 mm 20", ± 5 mm
Maximum Height	725 mm 29", ± 5 mm
Features	
Removable Head & Foot Panels	Available
Raised Backrest	Available
Raised Upper Leg	Available
DVT End Platform	Available (On Ratchet)
Four Section Metallic Perforated Tops	Available
125 mm High Grade Synthetic Single Wheel Castors, Two with Brake and Two without Brake	Available
Backrest Angle	0 - 75°, ± 2°
Knee Rest Angle	0 - 45°, ± 2°
Trendelenburg Angle	0 - 10°, ± 2°
Reverse Trendelenburg Angle	0 - 7°, ± 2°
Buffers at 4 Corners	Available
Heavy Duty I.V. Pole Provision at Four Corners	Available
Four separate screw mechanism for Backrest, Knee Rest, Trendelenburg & Reverse Trendelenburg Tilt & Height Adjustments	Available
Polymer Moulded Handle	Available
All M.S. Parts are Pre-treated and Powder Coated	YES
Polymer railing -(2 sets) (0360)	YES
Patient Load Bearing Capacity	100 kgs
Safe Working Load	135 kgs
Option of Head and Footboards	
Polymer Moulded Head and Footboards (0185) (Pair)	Available
Optional Accessories	
Four Section Mattress with 4" Thick PU Foam of 40 Density Covered with PVC (0321A) - 1900 mm x 865 mm, 75" x 34" (L x W)	YES
Oxygen Cylinder Cage (0166E)	YES
S.S. Heavy Duty I.V. Pole with Two Hooks (0164)	YES
S.S. Heavy Duty I.V. Pole with Four Hooks (0164C)	YES
Urine Bag Holder (0156)	YES
Moulded Chart Holder (0158A)	YES
S.S. Lifting Pole (0553)	YES
Orthopaedic Traction Pulley (0187A)	YES
Polymer railing single Set (Small) (0360A)	YES
Mattress arrester for comfort bed(0171C), 4 Nos	YES
Side Rails	Polymer Moulded Split Safety Rails (2 sets) 0360
Fixed polymer handle	With Fixed Polymer Moulded Handle (Set of Four) (in lieu of 0195 (Single Removable Poly. Mld. Handle))(0195B)
I.V. Rod with 2 hook	Telescopic IV pole with 2 hook. (0154B)
I.V. Rod with 4 hook	Telescopic IV pole with 4 hook. (0154D)
I.V. Rod with 2 hook	Utility I.V. Pole with 2 hook (0164B)
I.V. Rod with 4 hook	Utility I.V. Pole with 4 hook (0164D)

### Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070  
 Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com  
 Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



Cx4000 PLUS:01:25.10.17