

Cx4000B



Cx4000B - Comfort Bed (Fixed Height, Four Section & Two Way Longitudinal Tilt)

Comfort Bed having four sections with mild steel perforated sheet top. Bed with fixed height, having backrest, upper leg, trendelenburg and reverse trendelenburg positions on crank mechanism with polymer moulded handle. Stainless steel laminated head and foot boards and high grade synthetic body castors of 125 mm diameter, two with brake and two without brake.

Technical Data

Cx4000B

Overall Approximate Dimension	
Buffer to Buffer	2150 mm x 970 mm x 600 mm (L x W x H) 85" x 38" x 24"
Platform Dimension(**)	2060 mm x 920 mm (L x W) 81" x 36"
Mattress Platform	1950 mm x 790 mm (L x W) 77" x 31"
Buffer to Handle	2295 mm (L) 90"
Fixed Height	600 mm
Features	
Removable Head & Foot Panels	Available
Raised Backrest	Available
Raised Upper Leg	Available
Calf Elevation	Available
Metallic Perforated Tops	4 Section
125 mm High Grade Synthetic Castors, Two with Brake and Two without Brake	Available
Separate Screw Mechanism	Backrest, Upper Leg, Trendelenburg Tilt & Reverse Trendelenburg Tilt
Backrest Angle	Upto 65°
Knee Rest Angle	Upto 50°
Trendelenburg Angle	10°
Reverse Trendelenburg Angle	7°
Buffers at 4 Corners	Available
I.V. Pole Provision	Available
Polymer Moulded Handle * (0195)	Available
All M.S. Parts are Pre-treated and Epoxy Powder Coated	Yes
Patient Load Bearing Capacity	100 kgs
Safe Working Load	135 kgs
Option of Head and Footboards	
Polymer Moulded Head and Footboards (0183) (1 Pair)	Available
S.S. Laminated Head and Footboards (182G) (1 Pair)	Available
Optional Accessories	
Four Section Mattress with 4" Thick PU Foam of 40 Density Covered with PVC (0321) - 1980 mm x 865 mm (L x W) 78" x 34"	Yes
Oxygen Cylinder Cage (0166J)	Yes
S.S. Telescopic I.V. Pole (0154)	Yes
S.S. Telescopic I.V. Pole with 4 Hooks with S.S. Bracket (0154C)	Yes
S.S. Heavy Duty I.V. Pole (0164)	Yes
Heavy Duty I.V. Pole with 4 Hooks (0164C)	Yes
Urine Bag Holder (0156)	Yes
Moulded Chart Holder (0158A)	Yes
S.S. Lifting Pole (0553)	Yes
Orthopaedic Traction Pulley (0187A)	Yes
Full Length Collapsible Side Rails EPC (Pair) (0352)	Yes
Full Length Swing Down Type Side Rails EPC (Pair) (0354)	Yes
M.S. Collapsible Type Safety Side Rails Powder Coated (Pair) (0356A)	Yes
Collapsible Railings with Stainless Steel Tubes (0358) (Pair)	Yes

* Also available with Fixed Provision

** Also available in 1830 mm x 838 mm, 72" x 33"