

# Cx4000C

MIDMARK®

JANAK®



**Cx4000C - Comfort Bed** (Fixed Height, Two Section & Two Way Longitudinal Tilt)

Comfort Bed having two sections with mild steel perforated sheet top. Bed with fixed height having only backrest and trendelenburg / reverse trendelenburg position on crank mechanism with polymer moulded handle. Polymer moulded head and footboard and high grade synthetic body castors of 125 mm diameter, two with brake and two without brake.

## Technical Data

Cx4000C

Overall Approximate Dimension	
Buffer to Buffer	2050 mm x 970 mm x 600 mm (L x W x H) 81" x 38" x 24"
Platform Dimension(**)	2005 mm x 920 mm (L x W) 79" x 36"
Mattress Platform	1885 mm x 820 mm (L x W) 74" x 32"
Buffer to Handle	2275 mm (L) 90"
Fixed Height	600 mm
Features	
Removable Head & Foot Panels	Available
Raised Backrest	Available
Metallic Perforated Tops	2 Section
125 mm High Grade Synthetic Single Wheel Castors, Two with Brake and Two without Brake	Available
Separate Screw Mechanism	Backrest, Trendelenburg Tilt & Reverse Trendelenburg Tilt
Backrest Angle	Upto 65°
Trendelenburg Angle	10°
Reverse Trendelenburg Angle	7°
Buffers at 4 Corners	Available
I.V. Pole Provision	Available
Polymer Moulded Handle * (0195)	Available
All M.S. Parts are Pre-treated and Epoxy Powder Coated	Available
Patient Load Bearing Capacity	100 kgs
Safe Working Load	135 kgs
Option of Head and Footboards	
Polymer Moulded Head and Footboards (0183) (1 Pair)	Available
S.S. Laminated Head and Footboards (182G) (1 Pair)	Available
Optional Accessories	
Two Section Mattress with 4" Thick PU Foam of 40 Density Covered with PVC (0121) - 1920 mm x 865 mm (L x W) 75" x 34"	YES
Oxygen Cylinder Cage (0166J)	YES
S.S. Telescopic I.V. Pole with Two Hooks (0154)	YES
S.S. Telescopic I.V. Pole with Four Hooks with S.S. Bracket (0154C)	YES
S.S. Heavy Duty I.V. Pole (0164)	YES
Heavy Duty I.V. Pole with Four Hooks (0164C)	YES
Urine Bag Holder (0156)	YES
Moulded Chart Holder (0158A)	YES
S.S. Lifting Pole (0553)	YES
Orthopaedic Traction Pulley (0187A)	YES
Full Length Collapsible Side Rails (0352)	YES
Full Length Swing Down Type Side Rails EPC (Pair) (0354A)	YES
M.S. Collapsible Type Safety Side Rails Powder Coated (Pair) (0356A)	YES
Collapsible Railings with Stainless Steel Tubes (0358) (Pair)	YES

\* Also Available with Fixed Provision

\*\* Also available in 1830 mm x 838 mm, 72" x 33"