



1502A - Obstetric Labour Table Fully S.S.

- Overall approximate extended dimension: 1890 mm x 900 mm x 800 mm, 74" x 35" x 31.5" (L x W x H).
- Mattress platform dimension: 1800 mm x 655 mm, 71" x 26" (L x W).
- Three section S.S. top & S.S. tubular construction.
- Backrest adjustable on ratchet: 750 mm x 627 mm. 29.5" x 25" (L x W).
- Middle section with "U" cut: 350 mm x 655 mm. 14" x 26" (W x L).
- Leg section can be telescoped under the main frame completely for lithotomy positions: 595 mm x 700 mm. 23" x 27.5" (W x L).
- Foot operated gas spring assisted trendelenburg position.
- Complete with stainless steel head & side rails hand grips, a pair of rubber padded stainless steel lithotomy rods with rexine ankle straps and stainless steel bowl under 'U' cut. Size of railing: 195 mm x 685 mm. 8" x 27" (H x L).
- Supplied in SKD condition.

Optional Accessories

Lithotomy crutches in lieu of lithotomy rod (1552)

Rexine covered mattress of 25 mm (1") thickness (1551)

Stainless steel I.V. Rod - 2 hooks (0154) & 4 hooks (0154C)

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070

Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



1502A-1504A:02:25.10.17