



1512 - LDRP Bed

This motorized, versatile LDRP bed enables the obstetrician as well as the patient to perform their respective roles with comfort and ease. The bed is provided with handset for easy adjustments of backrest, height, trendelenburg and reverse trendelenburg.

Maternity & Gynaec Solutions



Minimum Height



Maximum Height



Pelvic Tilt



Trendelenburg



Reverse Trendelenburg

| Overall Approximate Dimension | |
|---|--|
| Buffer to Buffer | 2155 mm x 1010 mm (L x W) (85" x 40") (Bed when fully extended) |
| Frame Dimension | 2070 mm x 910 mm (L x W) (81" x 36") |
| Mattress Platform | 1800 mm x 650 mm (L x W) (71" x 26") |
| Minimum Height | 515 mm (20.2") |
| Maximum Height | 880 mm (35") |
| Features (Inbuilt) | |
| Patient Hand Control | Available |
| Central & Directional Locking | Available |
| Backrest Angle | 0 - 65° |
| Trendelenburg / Reverse Trendelenburg Angle | 0 - 25° |
| Buffers | Available |
| Telescopic I.V. Pole (0154) | Available |
| Battery Backup | Available |
| Polymer Moulded Head & Foot Bows | Available |
| Polymer Moulded Safety Side Rails (1 pair) | Available |
| Aluminium Lithotomy Crutches Mounted on Stainless Steel Rods (1 pair) | Available |
| Removable Stainless Steel Tray | Available |
| Three Section Mattress | Available |
| Dust Protective Base Cover | Available |
| 125 mm Diameter High Grade Synthetic Castors | Available |
| Mild Steel Parts Pre-treated and Powder Coated | Available |
| Safe Working Load | 180 kgs |
| Patient Load Bearing Capacity | 140 kgs |
| Electrical Specifications | |
| Power In (Mains Power Supply) | Nominal 230 V Ac 90 V Ac to 300 V Ac operating range with battery |
| Protection Voltage | 45 Hz to 55 Hz |
| Electrical Shock Protection | Class I |
| Degree of Shock Protection | Type B |
| Liquid Ingress Protection | IP 54 |
| Rechargeable Batteries | 2 x 12 Volt Sealed Lead Acid |
| Duty Cycle | 10% (Two minutes for every eighteen minutes) |

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



1512:02:25.10.17