

OT Solutions



1601 - Revolving Stool

- Tubular tripod based with S.S. revolving top.
- Height adjustable from approximately 465 mm (18") to 680 mm (27") by accurately machined screw mechanism.
- Legs fitted with rubber shoes.
- M.S. tubular construction.
- Pre-treated and powder coated finish.

1602 - Revolving Stool Fully Stainless Steel

- Tubular tripod based with stainless steel revolving top.
- Height adjustable from approximately 465 mm (18") to 680 mm (27") by accurately machined screw mechanism.
- Legs fitted with rubber shoes.
- S. S. tubular construction.

1811 - Two Tier Bowl Stand

- Approximate height: 850 mm. (33")
- Mild steel tubular frame mounted on three swivel, synthetic body castors 50 mm (2") wheel diameter without brake.
- Two stainless steel bowls of approximately 375 mm (15") diameter.
- Pre-treated and powder coated finish.