



1605A



1606A

1605A - Revolving Stool

- Tubular quadrapod based with S.S. revolving top with fixed cushion of 65 mm. (2.5") thickness
- Height adjustable from approximately 540 mm (21") to 755 mm (30") with cushion by accurately machined screw mechanism.
- Legs fitted with PVC shoes.
- M.S. tubular construction.
- Pre-treated and powder coated finish.

1606A - Revolving Stool Fully Stainless Steel

- Tubular quadrapod based with S.S. revolving top with fixed cushion of 65 mm. (2.5") thickness
- Height adjustable from approximately 540 mm (21") to 755 mm (30") with cushion by accurately machined screw mechanism.
- Legs fitted with PVC shoes.
- S. S. tubular construction.