



1605- Revolving Stool

- Tubular quadrapod based with S.S. revolving top
- Height adjustable from approximately 465 mm (18") to 680 mm (27") by accurately machined screw mechanism.
- Legs fitted with rubber shoes.
- M.S. tubular construction.
- Pre-treated and powder coated finish.
- S.S. Ring for footrest

1606- Revolving Stool Fully Stainless Steel

- Tubular quadrapod based with Stainless steel revolving top
- Height adjustable from approximately 465 mm (18") to 680 mm (27") by accurately machined screw mechanism.
- Legs fitted with rubber shoes.
- S. S. tubular construction.
- S.S. Ring for footrest

