

# Transfer Solutions



**JANAK**<sup>®</sup>  
Because we care.



Ergonomic Handgrip



Comfortable Backrest



Tubeless Polyurethane Tyres for Durability, Firm Grip & Safety Brake with Lever



Lightweight Aluminium Propelling Ring for Enhanced Grip



Detachable Armrest with Cushion



Anti-Bacterial Fabric



Non Corrosive Smooth Finishing Metal Frame



Detachable and Adjustable Footrest with 90° Swing on either side for Easy Patient Egress and Ingress



Adjustable and Foldable Footplates



1039



Seat Width  
460 mm



Overall Width  
635 mm



Folded Width  
290 mm



Overall Height  
935 mm



Turning Radius  
1820 mm



Overall Length  
1120 mm



Rear Wheels  
610 mm



Safe Load Bearing  
Capacity  
130 Kgs



Safety Belt  
Provided

Engineering Variation  $\pm 10$  mm

## Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

Check the nearest **Midmark India Experience Centre** at [www.midmark.in/miec](http://www.midmark.in/miec)

We are interested in growing relationships











TS-1039-1039A:00:12.04.18

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

# Transfer Solutions












**JANAK**<sup>®</sup>  
Because we care.

-  Ergonomic Handgrip
-  Comfortable Backrest
-  Tubeless Polyurethane Tyres for Durability, Firm Grip & Safety Brake with Lever
-  Detachable Armrest with Cushion
-  Anti-Bacterial Fabric
-  Non Corrosive Smooth Finishing Metal Frame
-  Detachable and Adjustable Footrest with 90° Swing on either side for Easy Patient Egress and Ingress
-  Adjustable and Foldable Footplates



## 1039A

|   |  |  |   |   |  |
|---|--|--|---|---|--|
| <br>Seat Width<br>460 mm | <br>Overall Width<br>635 mm | <br>Overall Height<br>935 mm  | <br>Overall Length<br>1120 mm | <br>Safe Load Bearing<br>Capacity<br>130 Kgs | <br>Safety Belt<br>Provided |
|   | <br>Folded Width<br>290 mm  | <br>Turning Radius<br>1820 mm | <br>Rear Wheels<br>305 mm     |   |  |

Engineering Variation  $\pm 10$  mm

### Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

Check the nearest **Midmark India Experience Centre** at [www.midmark.in/miec](http://www.midmark.in/miec)

We are interested in growing relationships



TS-1039-1039A:00:12.04.18

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark 'Midmark India' and the logo. Reproduction of this catalogue in any form is liable for legal action.