



1710 - Emergency & Recovery Trolley

- Overall approximate dimension: 1900 mm x 716 mm, 75" x 28" (L x W).
- Stretcher dimension: 1830 mm x 555 mm, 72" x 22" (L x W).
- Mattress platform size - 1790 mm x 555 mm, 70" x 22" (L x W).
- Backrest raise by ratchet up to 65°
- Trendelenburg Angle: 19°
- Reverse Trendelenburg Angle: 6°
- Two section top with removable X-Ray permeable stretcher, made of pre-laminated board supported on tubular frame.
- Height adjusted by foot operated high quality hydraulic pump from 660 mm (26") to 910 mm (36").
- Gas spring assisted trendelenburg / reverse trendelenburg positions.
- Backrest raised on ratchet.
- Four 125 mm (5") diameter castor wheels with high grade synthetic body, two with brake and two without brake.
- Buffers at four corners.
- Synthetic rubber covered pushing handles.
- Sliding X-Ray cassette holder.
- Oxygen cylinder cage and storage tray.
- Stainless steel telescopic I.V. Rod along with provision at four locations and swing-away type stainless steel side rails.
- Safety belts provided.
- Pre-treated and powder coated finish.
- Supplied in Semi Knocked Down condition (SKD).
- Safe working load of 165 kgs and patient load bearing capacity of 135 kgs.

1710 A - Emergency & Recovery Trolley

- Overall approximate dimension: 1900 mm x 855 mm, 75" x 34" (L x W).
- Stretcher dimension: 1830 mm x 710 mm, 72" x 28" (L x W).
- Mattress platform size - 1790 mm x 710 mm, 70" x 28" (L x W).
- Rest of the specifications as per 1710

Optional Accessories

Rexine covered mattress of 50 mm, (2") thickness for 1710 (1752D) and for 1710 A (1752F)