



1021 - Steel Stool (Fully Powder Coated)

- Overall approximate dimension: 305 mm x 305 mm x 510 mm, 12" x 12" x 20" (L x W x H)
- Base area size: 410 mm x 410 mm, 16" x 16" (L x W)
- Framework of mild steel, pressed bent angles and CRCA steel sheet top
- Legs fitted with PVC angular shoes.
- Pre-treated and powder coated finish

1022 - Steel Stool (S.S. Top & Powder Coated Framework)

- Overall approximate dimension: 305 mm x 305 mm x 510 mm, 12" x 12" x 20" (L x W x H)
- S.S. top size: 308 mm x 308 mm, 12" x 12" (L x W)
- Base area size: 410 mm x 410 mm, 16" x 16" (L x W)
- Superimposed S.S. top on framework of mild steel, pressed bent angles and CRCA steel sheet top
- Legs fitted with PVC angular shoes.
- Pre-treated and powder coated framework

1023 - Steel Stool (Fully S.S.)

- Overall approximate dimension: 305 mm x 305 mm x 510 mm, 12" x 12" x 20" (L x W x H)
- Base area size: 410 mm x 410 mm, 16" x 16" (L x W)
- Superimposed S.S. top on framework of stainless steel, pressed bent angles and CRCA steel sheet top
- Legs fitted with PVC angular shoes.

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



WCS-F:02:25.10.17



1603 - Single Step Stool

- Overall size: 510 mm x 305 mm, 20" x 12" (L x W).
- Step height: 232, 9".
- Mild steel tubular frame
- Mild steel top fitted with aluminium tread flats
- Legs fitted with PVC shoes.
- Pre-treated and powder coated finish

1604 - Double Step Stool

- Overall size: 603 mm x 510 mm x 452 mm, 24" x 20" x 18" (L x W x H).
- Step height: First 230 mm, 9" (H) & Second 452 mm, 20" (H).
- Mild steel tubular frame
- Mild steel top fitted with aluminium tread flats
- Leg fitted with PVC shoes.
- Pre-treated and powder coated finish

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



WCS-F:02:25.10.17