



Ax1012A Wardcare Adjustable Bedside Table

- Overall dimension: 960 mm x 420 mm, 38" x 16.5" (L x W). Top Size: 810 mm x 355 mm, 32" x 14" (L x W)
- Two sectional laminated top with railing on fixed section. Raised edge on one side tiltable section.
- Mild steel tubular telescopic stem, with geared S.S. handle for height adjustment from approximately 715 mm to 1080 mm, 28" to 43"
- Mild steel rectangular tubular base frame mounted on four twin wheel non rusting castors, 50 mm, 2" diameter without brake
- Pre-treated and powder coated
- Supplied in KDC

Ax1011 Wardcare Overbed Table

- Overall approximate dimension: 1120 mm x 450 mm x 980 mm, 44" x 18" x 38.5" (L x W x H)
- Laminated top of approximate 760 mm x 450 mm, 30" x 18" (L x W)
- Mild steel tubular base frame mounted on four swivel synthetic body castors 50 mm, 2" diameter without brake
- Pre-treated and powder coated finish
- Supplied in KDC

Ax1014 Stainless Steel Folding Food Table

- Overall approximate dimension: 560 mm x 400 mm x 280 mm, 22" x 16" x 11" (L x W x H)
- Stainless steel tray, machine pressed on edges, mounted on folding legs of stainless steel

Ward Care Solutions



0155B – Crib with Attachment (on Clamp Provision)

- Overall Crib Size: 760 mm x 380 mm x 300 mm (30" x 15" x 12") (L x W x H)
- Swings on two tubular rods fixed towards the leg side board of the bed
- Pre-treated & powder coated

**0155 - Crib with Attachment for Bx1000 Series (Specifications as per 0155B)*

**0122B - 3" PU Foam with PVC Rexine mattress*

1040 – Bed Elevator

- Overall Size: 650 mm x 300 mm (25.6" x 12") (L x W)
- Height adjustable from 980 mm to 1435 mm (38.5" to 56.5")
- Tubular frame mounted on two swivel castors of 100 mm (4") diameter, without brake
- Geared handle provided for easy lifting of the bed
- Hooks to hold and lift the bed either at head side or leg side of bed
- Pre-treated & powder coated



1051 A - Bed Side Screen (3Fold)

- Overall Size: 1680 mm x 2450 mm (66" x 96") (L x W)
- Mild steel tubular construction in three sections
- Mounted on six swivel twin wheel non rusting castors, 50 mm (2") diameter
- Middle frame 1220 mm (48") wide
- Side spans of 610 mm (24") each frame with four hooks
- Supplied with hooks and springs with curtains
- Pre-treated & powder coated
- Supplied in KDC

1069 - Ward Care - Cylinder Trolley (Push Type)

- Suitable for 1320 litre gas cylinder
- Mild steel tubular framework mounted on 125 mm (5") diameter castors
- Trolley provided with locking provision and stainless steel base
- Pre-treated & powder coated



Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070

Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



WCS-B;03;25.10.17