



Bx1003 - Attendant Bed cum Chair

Specifications

Overall Approximate Dimension	
Dimension from end to end	1850 mm x 785 mm, 73" x 31" L x W
Platform dimension (without armrest)	1850 mm x 625 mm, 73" x 25" L x W
Fixed height	525 mm, 21" H (without upholstery)
Upholstery thickness	75 mm, 3" H

Features

Backrest angle	Atleast 0° - 65°
Manual adjustment for back section and leg section	Available
Manual chair position	Available
Two armrest pivoted with backrest on each side	Available
Tray at bottom	Available
Powder coating	All M.S. parts are pre-treated and powder coated
Safe working load	135 kgs
Patient load bearing capacity	135 kgs



1038 - Commode Chair with ABS Seat & Cushion

- Overall approximate dimension: 512mm x 551 mm x 946 mm, 20" x 22" x 37" (W x D x H)
- Seat dimension: 405 mm x 454 mm, 16" x 18" (L x W)
- Removable Cushion and S.S. bowl
- S.S. tubular construction
- 125 mm diameter synthetic body castors, two with brake and two without brake
- Supplied in KDC