

Bx2000 PLUS



1. Conventional epoxy coated head and foot boards
2. Two section mattress with 4" thick PU foam of 40 density covered with PVC
3. Telescopic I.V. Pole with 2 hooks fully SS
4. Full length swing down type epoxy coated side rails
5. Removable polymer moulded handle
6. Urine bag holder



Bx2000 Plus with standard and additional accessories

Standard Accessories	Additional Accessories
1. Removable polymer moulded handle	1. Two section mattress with 4" thick PU foam of 40 density covered with PVC
	2. Conventional epoxy coated head and foot boards
	3. Full length swing down type epoxy coated side rails
	4. Telescopic I.V. Pole with 2 hooks fully SS
	5. Urine bag holder

Specifications:

- Bx2000 Plus has facility to raise the back section
- This section is operated by a crank mechanism provided at the foot end of the bed
- The raising section allow the nursing staff to change positions for various activities and clinical positions with relative ease
- Pre-treated and powder coated finish
- Supplied in SKD condition

* Option of fixed polymer moulded handle available in lieu of removable polymer moulded handle.

Bx2000 PLUS



JANAK[®]
Because we care.

Technical Data	Code	Bx2000 PLUS
Overall Approximate Dimension		
Platform Dimension		2055 mm x 915 mm / 91" x 36" (L x W)
Mattress Platform		1920 mm x 865 mm / 76" x 34" (L x W)
Fixed Height		600 mm H (optional height of 450 mm)
Features		
Raised Backrest		Available
Backrest Angle		Upto 65 degrees
I.V. Pole Provision		Available
Polymer Moulded Handle (Fixed / Removable)		Available
H-Type Legs with Rubber Shoes	0236E	Available
Option of Head and Foot Boards		
Polymer Moulded Head and Foot Boards	0183	Available
S.S. Laminated Head and Foot Boards	0182	Available
Conventional Head and Foot Boards	0181	Available
MS Head & Foot Boards	0181E	Available
Optional Accessories		
Two Section Mattress (with 4" thick PU foam of 40 density covered with PVC)	0121	Yes
Two Section Mattress (with 3" thick PU foam of 32 density covered with PVC)	0121C	Yes
Mattress Arrestor (set of two, location needs to be specified)	0171	Yes
Telescopic I.V. Pole - Two Hooks	0154	Yes
Telescopic I.V. Pole - Four Hooks	0154C	Yes
Urine Bag Holder	0156	Yes
Moulded Chart Holder	0158A	Yes
H-Type Legs 910mm (36") Width and 600mm (18") Height with Rubber Shoes for Basic Bed (For Bx2000E,Bx4000E)		Yes

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships

