

Bx4000 Basic Fowler Bed



Basic Folwer Bed with conventional head and foot panels (0181)



Basic Folwer Bed with stainless steel laminated head and foot panels (0182G)



Basic Folwer Bed polymer moulded head and foot panels (0183)

Specifications:

- The Basic Folwer Bed has facility to raise the back and upper leg sections
- These sections are operated by crank mechanisms provided at the foot end of the bed
- The raising sections allow the nursing staff to change positions for various activities and clinical positions with relative ease
- Pre-treated and powder coated finish
- Supplied in SKD condition

* Option of fixed polymer moulded handles (set of two) available in lieu of single removable polymer moulded handle.

Bx4000 Basic Fowler Bed



Technical Data	Bx4000
Overall Approximate Dimension	
Bed Frame Dimensions (L x W)	2070 mm x 920 mm / 82' x 36", ± 10mm
Mattress Platform (L x W)	1980 mm x 790 mm / 78" x 34", ± 10mm
Buffer to Handle (L)	2275, 90", ± 10mm
Fixed Height with H-leg with rubber bush	600 mm H (optional height of 450 mm). ± 5mm
Fixed Height with H-leg with high grade 125 mm synthetic castor	600 mm H (optional height of 450 mm). ± 5mm
Buffer to Buffer (when castors provided) (L x W)	2140 mm x 970 mm , 84' x 39", ± 10mm
Features	
Raised Backrest, Upper Leg, Calf Elevation	Available
Backrest Angle	Upto 65°, ± 2°
Knee Rest Angle	Upto 40°, ± 2°
I.V. Pole Provision	Available
Polymer Moulded Handle (Fixed / Removable)	Available
H-Type Legs with Rubber Shoes (0236)	Available
Option of Head and Footboards	
Polymer Moulded Head and Foot Boards (0183)	Available
S.S. Laminated Head and Foot Boards (0182G)	Available
Conventional Head and Foot Boards (0181)	Available
Safe Working Load	135 kg
Patient Working Load	100 kg
Optional Accessories	
Four Section Mattress with 4" thick PU foam of 40 density covered with PVC (0321)	Yes
Mattress Arrestor set of two (0171)	Yes
Telescopic I.V. Pole - Two Hooks (0154), Four Hooks (0154C)	Yes
Urine Bag Holder (0156)	Yes
Mosquito Net Poles (0157D)	Yes
Moulded Chart Holder (0158A)	Yes
S.S. Lifting Pole (0553)	Yes
Orthopaedic Traction Pulley (0187A)	Yes
Full Length Swing Down Type Side Rails (Pair) (0354A)	Yes
M.S. Collapsible Type Safety Side Rails - Powder Coated (Pair) (0356A)	Yes
Collapsible Railing with Stainless Steel Tube (Pair) (0358)	Yes
Full length Stainless Steel Collapsible Side Rails (0352)	Yes
H-Type Legs 910mm, 36" (W), 450mm, 18" (H) with Rubber Shoes (0236A)	Yes
H-Type Legs 910mm, 36" (W) with 125mm High Grade Synthetic Body Castors, two with brake and two without brake (0236B) - for 600mm, 24" Height	Yes
H-Type Legs 910mm, 36" (W) with 125mm Synthetic Body Castors, two with brake and two without brake (0236C) - for 600mm, 24" Height	Yes
H-Type Legs 910mm, 36" (W) with 125mm High Grade Synthetic Body Castors, two with brake and two without brake (0236D) - for 450mm, 18" Height	Yes
H-Type Legs 910mm, 36" (W) with 125mm Synthetic Body Castors, two with brake and two without brake (0236E) - for 450mm, 18" Height	Yes
Crib with Attachment with Clamp (0155B) with 3" PU foam mattress with PVC rexine mattress	Yes
Rubber Buffers at Corners (when castors provided)	Yes

** While attaching the Crib, user has to remove the Polymer Moulded Bow (0183) at foot end to use it.
(For Bx4000 the upper and lower leg to be flat when using the Crib)

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070

Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



Bx4000;04;17.07.18