



0104H

Technical Specifications

Buffer to Handle*	1720 mm, 68"
Buffer to Buffer*	1500 mm x 820 mm, 59" x 32" (L x W)
Platform Dimension*	1450 mm x 775 mm, 57" x 31" (L x W)
Mattress Platform*	1335 mm x 665 mm, 53" x 26" (L x W)
Height*	480 mm to 720 mm, 19" to 28"
Backrest **	0° - 70°
Trendelenburg **	0° - 10°
Reverse Trendelenburg **	0° - 8°
Safe Working Load	110 kgs
Patient Load Bearing Capacity^	90 kgs

Standard Features

- Four corner safety rubber buffer of 100 mm diameter
- SS tubular, telescopic, collapsible side rails
- High grade non floor-staining synthetic material castor with integrated thread guards.
- SS laminated head and foot boards
- Telescopic IV pole with 2 hooks
- Removable handle
- All MS parts epoxy powder coated

Add On Accessories

- Two section mattress with 3" thick PU foam of 32 density covered with PVC
- Chart Holder
- Urine Bag Holder
- Oxygen Cylinder Cage with Epoxy Powder Coated
- SS telescopic IV pole with 2 hooks
- SS heavy duty IV pole with 2 hooks

Safety Standards

General Safety for medical beds
IEC 60601-2-52

*± 10 mm Engineering Variation and ** ± 2° Engineering Variation ^ When Bed is configured with standard accessories

 *± 10 mm Engineering Variation

We are interested in growing relationship



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

0104H:05:08.01.24