

# ICU Solutions



**JANAK®**  
Because we care.



## 9003 Crash Cart

Overall approximate dimension:  
960 mm x 500 mm x 1635 mm, 38" x 20" x 64" (L x W x H)\*

SS tubular frame and SS shelves

Two rows with six multicoloured removable bins and polystyrene lockable storage units with three drawers having containers of different sizes

125 mm high grade synthetic castors two with brakes and two without brakes

Complete with powder coated oxygen cylinder holder, catheter holder, cardiac massage board, SS IV pole

Safe working load: 50 kgs

Supplied in KDC condition



## 9001 Drug Trolley

Overall approximate dimension:  
680 mm x 410 mm x 1570 mm, 27" x 16" x 62" (L x W x H)\*

Large bin dimension:  
250 mm x 165 mm x 125 mm, 10" x 6" x 5" (L x W x H)\*

Small bin dimension:  
150 mm x 110 mm x 75 mm, 6" x 4" x 3" (L x W x H)\*

30 multicoloured removable bins including 12 nos. of large and 18 nos. of small bins

SS shelf at the top and MS tray at the bottom for medicinal storage

125 mm high grade synthetic castors two with brakes and two without brakes

Safe working load: 50 kgs

Supplied in KDC condition

\*± 10 mm Engineering Variation

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.



# ICU Solutions



## 9003A Crash Cart

Overall approximate dimension:  
1090 mm x 400 mm x 1650 mm, 43" x 16" x 65" (L x W x H)\*

MS tubular frame and MS shelves on top and bottom

Most frequently used SS centre shelf

Two rows with six multicoloured removable bins and polystyrene lockable storage units with three drawers having containers of different sizes

125 mm high grade synthetic castors two with brakes and two without brakes

Complete with powder coated oxygen cylinder holder, cardiac massage board, SS IV pole

Safe working load: 50 kgs

Supplied in KDC condition

\*± 10 mm Engineering Variation



### Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



9003-9001-9003A:02:20.07.21