

Maternal & Gynaec Solutions

1512 LDRP Bed



This motorized, versatile bed for Labour, Delivery, Recovery, Postpartum enables the obstetrician to deliver the promise of safe birthing and recovery with a carefree experience. The bed provides comfort and ease to the patient as well as the obstetrician to assist in the hospital delivery. The bed is provided with wired handset control for easy adjustment of backrest, height, trendelenburg and reverse trendelenburg functions.



We are interested in growing relationship





Minimum Height



Maximum Height



Pelvic Tilt



Trendelenburg



Reverse Trendelenburg

Technical Specifications

Overall Dimension (Buffer to Buffer)*	2160 x 1010 mm (LxW)
Frame Dimension*	2070 x 910 mm (LxW)
Mattress Platform Dimension*	1805 x 650 mm (LxW)
Backrest**	0° - 65°
Trendelenburg & Reverse Trendelenburg**	0° - 25° & 0° - 23°
Height*	525 mm to 892 mm
Safe Working Load	180 Kgs
Patient Load Bearing Capacity^	140 Kgs

Standard Features

- Wired handset control for patients
- Telescopic IV pole with 2 hooks
- Three section mattress with pink colour laminate
- Aluminium lithotomy crutches mounted on stainless steel rods (1 pair)
- Dust protective base cover
- Polymer moulded safety side rails
- Polymer moulded head and foot boards
- 125 mm high grade synthetic twin castors with central and directional locking mechanism
- Four corner safety rubber buffer of 125 mm diameter
- Removable SS tray
- All MS parts pretreated and epoxy powder coated
- Battery backup

Safety Standards


• General safety for medical electrical equipment	IEC 60601-1-1
• EMI/EMC	IEC 60601-1-2
• General safety for medical bed	IEC 60601-2-52

Electrical Specifications


Power In (Operating Range)	100 VAC to 240 VAC 50/60 Hz, Max 5A
Liquid Ingress Protection	IPX6
Electrical Shock Protection	Class 1
Degree of Shock Protection	Type B
Duty Cycle	10% (two minutes for every eighteen minutes)

* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When bed is configured with standard accessories. #Factory fitted.

Add-On Accessories



Accessories Rack



Accessories rack with installed accessories

- Telescopic IV pole with 4 hooks

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1512:06:06.06.23



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark" India, and the logo. Reproduction of this catalogue in any form is liable for legal action.