

Transfer solutions



JANAK[®]
Because we care.



1036

1036 – Stainless steel Rigid Wheelchair with Cushioned Seat and Back

- Overall approximate dimension:
1120 mm x 670 mm W x 920 mm H 44" x 26" x 36" (D x W x H)*
- S.S. Tubular framework fitted with cushioned seat and back
- Adjustable aluminium footrest
- Two solid rubber tyred bicycle wheels with brake and self-propelling S.S. hoops
- Two swivel castors in front of 200 mm diameter (8")*
- S.S. tubular construction
- Safety belt provided
- Moulded chart holder provided
- Weight bearing capacity: 135 Kgs
- Upholstery colour: Pista



1033A

1033A – Wheelchair

- Overall approximate dimension:
1120 mm x 670 mm W x 920 mm H 44" x 26" x 36" (D x W x H)*
- Mild steel tubular framework fitted with cushioned seat and back
- Moulded chart holder provided
- Adjustable aluminium footrests
- Two solid rubber tyred bicycle wheels with brake and self-propelling S.S. hoops
- Two swivel castors, 200 mm* diameter in front
- Safety belt provided
- Pre-treated and powder coated
- Supplied in Knock Down Condition (KDC)
- Weight bearing capacity: 135 Kgs
- Upholstery colour: Pista

*(± 10 mm Engineering Variation)



Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



1033A-1036:04:03.06.19