



## 1807 – Health Check Examination Table with Upholstered Top

- Overall approximate dimension: 1835 mm x 610 mm x 826 mm, 72" x 24" x 32.5" (L x W x H).
- Mild steel tubular construction, machine pressed double bend mild steel sheet.
- Two section top with fixed upholstery.
- Headrest adjusted on ratchet with four elevation angles.
- Legs fixed with rubber shoes.
- Pre-treated and epoxy powder coated.
- Supplied in KDC.
- Safe working load 140 kgs and patient load bearing capacity of 135 kgs.
- Upholstery Colour: Pista