

Exam Room Solutions

1803H

Manual Examination Table



1803H

Designed with you in mind.

The 1803H manual examination table is the outcome of research focusing on design, comfort, and storage. Featuring six drawers and a cabinet, it provides efficient storage at the point of care for items such as gowns and drapes, along with well-organized space for smaller supplies and bulk items. This configuration promotes flexibility in room arrangement and workflow, facilitating easy access to essential supplies.

Stylish. Comfortable. Easy to clean.

The 1803H upholstered top boasts a seamless design that not only feels as soft and comfortable as it appears but also minimizes areas where dirt and bacteria can gather. This seamless upholstery design simplifies cleaning, making it quicker, more efficient, and more effective.

We are interested in growing relationship





Thoughtful Design

Oval design to differentiate your set up



Unmatchable Comfort

As soft and comfortable as it looks



Adequate Storage Space

Providing easy access to essential supplies



Patient Positioning

For better access during examination

Optional Accessories



Sliding Step Stool

Aluminum grip on the step stool for easy ingress/egress



Paper Roll Holder

To minimize surface contamination



Adjustable BP Apparatus Tray

Height and multiple position adjustable BP apparatus tray

Specifications

Safe Working Load: 170 kg

Patient Working Load: 135 kg

Overall Dimensions: 1910 mm x 720 mm x 785 mm L x W x H

Mattress Thickness: 65 mm

Back Support Angle: Minimum = 0°; Maximum = 70°

Backrest dimensions: 560 mm x 670 mm LxW

Metal Structure Length: 1230 mm

Metal Structure Width: 550 mm

Dimensions of Drawer under gas spring: 410 mm x 335 mm x 80 mm LxWxH

Dimensions of other Drawers: 410 mm x 335 mm x 140 mm LxWxH

Drawer Load Capacity: 2 kg

Cabinet Dimensions: 368 mm x 365 mm x 515 mm WxHxD

Cabinet Load Capacity: 5 kg

Sliding stool usable area dimensions: 285 mm x 360 mm LxW

Max Width of the paper roll holder: 400 mm

CE *± 10 mm Engineering variation. ** ± 5° Engineering variation.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th Floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070.
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1901:03:11.09.24



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.