



AX2001 Supreme Care Bedside Locker

- Overall size (with railing): 495 mm x 475 mm x 865 mm, 19.5" x 18.7" x 34" (LxWxH)
- It is constructed from superior compact laminate, reinforced with a robust aluminum profile with zinc anodizing for extended durability
- Both the drawer front and cabinet door are crafted from laminated materials and overlaid with PVC foil that is resistant to scratches and UV rays, with a thickness of 400 microns
- The drawer measures 355 mm x 380 mm x 90 mm (LxWxH) and is equipped with smooth slides
- The dimensions of the open storage space is 380 mm x 380 mm x 135 mm
- The dimensions of the cabinet is 380 mm x 380 mm x 450 mm
- The base of the drawer features dual-wheel castors, each with a diameter of 50 mm, made from premium polyurethane
- The half-railing design provides additional support



AX1010 Super Deluxe Bedside Locker

- Overall approximate size: 400 mm x 400 mm x 820 mm (L x W x H)* (without file holder and towel holder)
- Overall approximate size: 490 mm x 450 mm x 820 mm (L x W x H)* (with file holder and towel holder)
- Drawer size: 150 mm x 350 mm x 365 mm (H x W x D)* approx. fitted with very smooth sides
- Front top drawer cover and cabinet door cover is made from ABS and are laminated with scratch resistant and UV resistant 350 micron PVC sheet
- Body is made from machine pressed one piece CRCA MS sheet which is pre-treated and with powder coated finish
- It is equipped with file holder and towel holder
- Doors and drawer are provided with epoxy coated aluminium handle
- Base is fitted with four swivel castors of 50 mm diameter made from synthetic material without brake

AX1010A Super Deluxe Bedside Locker

- Lock provided on drawer
- Rest of the specifications are the same as AX1010



AX1009

AX1009 Deluxe Bedside Locker

- Overall approximate size: 400 mm x 400 mm x 820 mm, 16" x 16" x 32" (L x W x H)*
- Overall approximate size from handle to buffer: 454 mm, 18" (D)*
- Drawer size: 90 mm x 355 mm x 380 mm, 3.5" x 14" x 15" (H x W x D)*
- Locker made from machine pressed CRCA MS Sheets enclosed on three sides which is pre-treated and powder coated finish
- Having one box and one drawer, fitted with membrane top with raised edges, drawer and cabinet door on front side and buffers at rear to avoid collision
- Base is fitted with four swivel castors of 50mm (2") diameter made from synthetic material without brake

AX1009B Deluxe Bedside Locker

- Lock provided on drawer
- Rest of the specifications are the same as AX1009



AX1006

AX1006 Wardcare Bedside Locker with SS Top

- Overall size: 410 mm x 420 mm x 820 mm, 16" x 17" x 32" (L x W x H)*
- Drawer size: 100 mm x 350 mm x 390 mm, 4" x 14" x 15" (H x W x D)*
- Locker made from machine pressed CRCA steel sheets enclosed on three sides having one box and one drawer with side
- Fitted with superimposed stainless steel top with raised edge on three sides
- Mounted on four swivel, synthetic body castors, 50 mm, (2") wheel diameter without brake
- Buffers at the rear
- Pre-treated and powder coated finish

AX1006A Wardcare Bedside Locker

- Lock provided on drawer
- Rest of the specifications are the same as AX1006



AX1001

AX1001 Wardcare Bedside Locker with SS Top

- Overall size: 410 mm x 410 mm x 835 mm, 16" x 16" x 33" (L x W x H)*
- Cabinet box of size 312 mm x 335 mm x 330 mm, 12" x 13" x 15" (H x W x D)*
- Machine pressed CRCA steel sheet locker box with door
- Fitted with stainless steel top with raised edges on three sides
- Mild steel tubular legs fitted with PVC shoes
- Pre-treated and powder coated finish
- Supplied in KDC

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

CE *± 10 mm Engineering Variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

WCS-A:10:21.05.24

